

Product Sheet: Anti-Hepatitis B virus core antigen IgG fraction (monoclonal)

Catalog Number	HBM-021-5
Product Id	HBM-021-5
Product Name	Anti-Hepatitis B virus core antigen IgG fraction (monoclonal)
Temperature	-20°C
Unit	100µg
Price	295
Availability	Available
Category	Infectious Agents
Subcategory	Hepatitis B Virus
Description	<p>Anti-hepatitis B core antigen IgG fraction is a mouse monoclonal antibody (Clone 2C11/2) to native HBcAg particles.</p> <p>The antigen used for immunization is highly purified HBcAg produced by recombinant E. coli cells.</p> <p>It has been purified using Protein G.</p>
Isotype	IgG-2a
Mol Weight	N/A
Purity	N/A
Storage	Store at -20°C.
Stability	At least 1 year at -20°C.
Biological Activity	The antibody is reactive in ELISA and also immunoprecipitation of native hepatitis B core particles.
Formulation	Solution at 0.50 mg/ml in 0.1M Tris-HCl (pH 7.5)
References	N/A

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
 184 Park Drive, Milton Park
 Abingdon, UK
 T: +44 (0)1235 828 200
 F: +44 (0) 1235 820 482

 **North America**
 1035 Cambridge Street,
 Cambridge, MA 02141
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

 **Germany**
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880

 **Switzerland**
 Centro Nord-Sud 2E
 CH-6934 Bioggio-Lugano
 T: +41(0) 91 604 55 22
 F: +41(0) 91 605 17 85