

Product Sheet: Anti-Hepatitis C virus NS5 antigen IgG fraction (monoclonal)

Catalog # HCM-131-5

Temperature	-20°C
Unit	100µg
Category	Infectious Agents
Subcategory	Hepatitis C Virus
Description	Mouse monoclonal antibody recognizing the NS5 antigen (Clone 2F6/G11) spanning N2054-C2295 of the HCV (subtype 1a) sequence. It has been produced using highly purified NS-5 antigen.
Isotype	IgG-1
Mol Weight	N/A
Purity	N/A
Storage	Store at -20°C.
Stability	At least 1 year at -20°C.
Biological Activity	This antibody is reactive in ELISA and dot blot tests. IgG Subtype: IgG-1.
Formulation	Solution at 0.50 mg/ml in 0.1M Tris-HCl (pH 7.5)
References	N/A

All products are for research use ONLY in laboratory animals and in vitro testing. NOT for diagnostic or human use.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85