

HERM-1811-5

Anti-human Endogenous Retrovirus type K (HERV K) Envelope protein. IgG fraction (monoclonal)

Product Id	HERM-1811-5
Product Name	Anti-human Endogenous Retrovirus type K (HERV K) envelope protein. IgG fraction (monoclonal)
Temperature	-20 °C
Unit	100 µg
Description	IgG fraction from mouse monoclonal antibody (clone 215/H4) obtained by immunization of mice with HERV K envelope antigen (79.5 kDa). The IgG fraction was purified using Protein G-Sepharose.
Isotype	IgG2a
Storage	Store at -80 °C.
Stability	Stable for at least one year. Avoid repeated freezing and thawing.
Biological Activity	This antibody can be used in ELISA assays (1:8,000 dilution).
Formulation	Solution at 1.0 mg/ml in PBS.
References	<ol style="list-style-type: none"> Toenjes, R.R. et. al. (1997) Virology 233, 280 - 291. Barbulescu, M. et. al. (1999) Curr. Biol. 9, 861 - 868.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85