

EZH2 (Gene ID 2146) qPCR Primer Pairs

Specifications

Catalog No. **Description**

HP207764 qSTAR qPCR primer pairs against Homo sapiens gene EZH2

Forward Sequence: GACCTCTGTCTTACTTGTGGAGC

Reverse Sequence: CGTCAGATGGTGCCAGCAATAG

Accession No. NM_001203247, NM_001203248, NM_001203249, NM_004456, NM_152998

Synonyms ENX-1; ENX1; EZH1; EZH2b; KMT6; KMT6A; WVS; WVS2

Component 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)

Quality Control The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage & Stability The primer mix is stable for one year from date of shipping when stored at -20 degrees C

* Primer pairs are not supplied with RT-PCR reagents

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85