

PDCD1 (Gene ID 5133) qPCR Primer Pairs

Specifications

Catalog No.	Description
HP208236	qSTAR qPCR primer pairs against Homo sapiens gene PD-1 / PDCD1

Forward Sequence: AAGGCGCAGATCAAAGAGAGCC

Reverse Sequence: CAACCACCAGGGTTTGGAAGT

Accession No. NM_005018

Synonyms CD279; hPD-1; hPD-I; hSLE1; PD-1; PD1; SLEB2

Component 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)

Quality Control The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage & Stability The primer mix is stable for one year from date of shipping when stored at -20 degrees C

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85