

MAP1LC3B (Gene ID 81631) qPCR Primer Pairs

Specifications

Catalog No.	Description
HP214570	qSTAR qPCR primer pairs against Homo sapiens gene MAP1LC3B

Forward Sequence: GAGAAGCAGCTTCCTGTTCTGG

Reverse Sequence: GTGTCCGTTACCAACAGGAAG

Accession No. NM_022818

Synonyms ATG8F; LC3B; MAP1A/1BLC3; MAP1LC3B-a

Component 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)

Quality Control The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage & Stability The primer mix is stable for one year from date of shipping when stored at -20 degrees C

UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

North America

23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com
AMS Biotechnology