

-20°C**ANTI-H. PYLORI UREASE (SMALL SUBUNIT) IgG FRACTION
(MONOCLONAL)**Catalog number: **HPM-5021-5 (100µg)** Lot #:

Description and Biological Activity: Mouse monoclonal antibody raised from highly purified H. pylori urease small subunit (Clone 2E4/F5). It has been purified using protein G.

Isotype: IgG-1.

Biological Activity: This antibody may be used on an ELISA or Western blot analysis (1/100-1/2500 dilution).

Storage: Store at -20°C.

Stability: Stable at least one year at -20°C.

Formulation: Solution at 0.50 mg/ml in 0.1M Tris-HCl (pH 7.4). Some of the liquid in the vial could have evaporated with changes in the final volume. However, the mass of the protein is still inside the vial.

For Research ONLY. Not for Human or Diagnostic Use.