

Catalog Number	HTIM-981-5
Product Id	HTIM-981-5
Product Name	Anti-HTLV-I gag p19 antigen IgG fraction (monoclonal)
Temperature	-20°C
Unit	100µg
Price	425
Availability	Available
Category	Infectious Agents
Subcategory	Human T-Lymphotropic Virus I
Description	IgG fraction from mouse monoclonal antibody obtained by immunization with highly purified HTLV-I gag p19 antigen (Clone 6E10/B8)
Isotype	IgG-2b
Mol Weight	N/A
Purity	N/A
Storage	Store at -20°C.
Stability	At least 1 year at -20°C.
Biological Activity	Reacts in ELISA and Western blot (1/100-1/1000 dilution).
Formulation	Solution at 0.10 mg/ml in 0.1M Tris-HCl (pH 7.5)
References	N/A

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, UK
 T: +44 (0)1235 828 200
 F: +44 (0) 1235 820 482

 **North America**
 1035 Cambridge Street,
 Cambridge, MA 02141
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

 **Germany**
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880

 **Switzerland**
 Centro Nord-Sud 2E
 CH-6934 Bioggio-Lugano
 T: +41(0) 91 604 55 22
 F: +41(0) 91 605 17 85