

## Human adult normal paraffin tissue sections

**Quantity:** 5 slides

**Storage Conditions:** Cool dry at room temperature. Keep at 4°C for long term storage. One year from the date of receipt under proper storage condition.

**Applications:** Human adult normal paraffin tissue sections are ideal for:

- Cellular localization of specific mRNA and protein expression in human normal tissues.

**Quality Control:** At least one of the tissue slides from each lot was stained with H & E to ensure the quality.

**Description:** A single human adult normal tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections are 5 µm in thickness and mounted on each positively charged glass slides.

**Shipping:** Room temperature.

**Intended use of the products:** All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. AMSBIO products may not resell, modified for resale or used to manufacture commercial products without written approval from AMSBIO.

**Limited product warranty:** AMSBIO has a total quality assurance policy. If any product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits AMSBIO liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not be liable for any incidental, direct, special or consequential damages, even if it is aware of the possibility of such damages.