

3X NASBA Reaction Buffer

Catalog Number: NEC-B-24


Storage Temperature: -20°C

Unit: 200µl (24 Reactions)

Description: This reaction buffer is to be used in conjunction with the liquid NASBA Enzyme Cocktail (NEC-1-24) in NASBA reactions. It requires the addition of nucleotides by the user.

Composition Per Reaction:

- 40 mM TrisHCl, pH 8.5 @ 25°C
- 12 mM MgCl₂
- 70 mM KCl
- 10 mM DTT
- 15% Dimethyl Sulfoxide



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

North America

23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com
AMS Biotechnology