

Product Sheet: Rabbit anti-Trout Growth Hormone II IgG fraction (polyclonal)

Catalog Number	PA-223-9
Product Name	Rabbit anti-Trout Growth Hormone II IgG fraction (polyclonal)
Temperature	-20°C
Unit	1mg
Category	Growth Factors, Cytokines, Receptors and their Antibodies
Subcategory	Polyclonal Antibodies (Growth Factors/Cytokines/Hormones)
Description	Anti-Trout growth hormone II IgG fraction is a polyclonal antibody raised in rabbit against purified Trout Growth Hormone II. It is purified using Protein G.
Isotype	N/A
Mol Weight	N/A
Purity	N/A
Storage	Store at -20°C.
Stability	At least 1 year at -20°C.
Biological Activity	This polyclonal antibody recognizes native Trout Growth Hormone II on an ELISA format assay.
Formulation	IgG fraction. Frozen liquid. 1.0 mg/ml. Formulated in 0.1 M Tris-HCl (pH 7.5)
References	N/A

All products are for research use ONLY in laboratory animals and in vitro testing. NOT for diagnostic or human use.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 617 945 5033 or
T: +1 800 987 0985
F: +1 617 945 8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85