



Cat# **RS-AAV1-ET**

AAV1 Reference Standards (Empty Capsids)

Type of Standards: AAV Full Capsids

Serotype: AAV1

Volume: 100 µl

Serotype: AAV1

Storage Buffer: PBS with 0.01% Pluronic F-68

Storage Conditions: -80°C

Shipping: Dry ice

Viral Genome Titer: 2.26×10^{12} VP/mL

Percentage Of Empty Capsids: 92.60%

For Research Use Only. Not for Use in Therapeutic or Diagnostic Procedures.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85