



Cat# **RS-AAV1-FL**

### **AAV1 Reference Standards (Full Capsids)**

Type of Standards: AAV Full Capsids

Serotype: AAV1

Promoter: CMV

Gene: GFP

Volume: 100 µl

Serotype: AAV1

Storage Buffer: PBS with 0.01% Pluronic F-68

Storage Conditions: -80°C

Shipping: Dry ice

Viral Genome Titer:  $5.33 \times 10^{11}$  GC/mL

Percentage Of Full Capsids: 61.00%

For Research Use Only. Not for Use in Therapeutic or Diagnostic Procedures.

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

 **UK & Rest of the World**  
184 Park Drive, Milton Park  
Abingdon OX14 4SE, UK  
T: +44 (0)1235 828 200  
F: +44 (0) 1235 820 482

 **North America**  
1035 Cambridge Street,  
Cambridge, MA 02141  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

 **Germany**  
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880

 **Switzerland**  
Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
T: +41(0) 91 604 55 22  
F: +41(0) 91 605 17 85