



Cat# **RS-AAV2-ET**

AAV2 Reference Standards (Empty Capsids)

Type of Standards: AAV Full Capsids

Serotype: AAV2

Volume: 100 µl

Serotype: AAV2

Storage Buffer: PBS with 0.01% Pluronic F-68

Storage Conditions: -80°C

Shipping: Dry ice

Viral Genome Titer: 1.27×10^{12} VP/mL

Percentage Of Empty Capsids: 99.50%

For Research Use Only. Not for Use in Therapeutic or Diagnostic Procedures.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85