

Cat # <small>Ambsbio</small>	RS-AAV5-FL
Description	AAV5 Reference Standards (Full Capsids)
Volume	100 µl
Promoter	CMV
Serotype	AAV5
Storage Buffer	PBS with 0.01% Pluronic F-68
Storage Conditions	-80°C
Shipping	Dry ice
Viral Particle Titer	2.60x10 ¹¹ GC/mL
Percentage of Full Capsids	67.30%

Use Restrictions

For Research Use Only. Not for Use in
Therapeutic or Diagnostic Procedures.

Storage Conditions



-80°C

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

North America 
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

Europe 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0)72 8080244
F: +31 (0)72 8080142

Switzerland 
Via Lisano 3,
(CP.683)
CH-6900
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

Germany 
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880