

CERTIFICATE OF ANALYSIS

LOT COMPOSITION:

- ▶ **Number of donor(s) in lot** 19
- ▶ **Gender of donor(s):** M/F

CATALOG #: SER-PM**Lot #:** 032911**Protein Concentration:** 15mg/ml**DESCRIPTION:** Human Platelet Membrane Cell Lysate

Human Platelet Membranes were derived from platelet rich plasma of pooled adult donors. All donors are fully consented from consented adult donors. The platelet membrane protein is isolated via homogenization and ultracentrifugation. The product is shipped frozen on dry ice. Store at -20°C or greater. (-80°C is preferred for long-term storage)

QUALITY CONTROL

Each lot is tested via PCR and found non-reactive to viral DNA from HIV-1, HIV-2, HTLV I, HTLV II, hepatitis B, hepatitis C and syphilis. However, no known test can offer complete assurance that the viruses that cause HIV-1, HIV-2, HTLV I, HTLV II, hepatitis B, hepatitis C and syphilis are not present. Since we cannot test all the pathogens, please treat the culture as a potential infectious reagent at Biosafety Level 2 as recommended by the CDC/NIH Manual.

These results have been verified.

SIGNATURE ON FILE

08/22/2013

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85