

<b>Name:</b> HABP2 mouse monoclonal antibody,clone OTI8D6 <b>Product Data Sheet - TRUEMAB</b>	<b>Catalog:</b> TA809299
--	--------------------------


<b>Components:</b>	• HABP2 mouse monoclonal antibody,clone OTI8D6 (TA809299)
<b>Amount:</b>	100ul
<b>Immunogen:</b>	Human recombinant protein fragment corresponding to amino acids 33-330 of human HABP2 (NP_004123) produced in E.coli.
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG2a
<b>Species Reactivity:</b>	Human
<b>Guaranteed Applications:</b>	WB, IHC
<b>Suggested Dilutions:</b>	WB 1:500, IHC 1:2000,
<b>Concentration:</b>	1.00mg/ml
<b>Buffer:</b>	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Purification:</b>	Purified from cell culture supernatant by affinity chromatography
<b>Storage Condition:</b>	Shipped at -20C or with ice packs. Upon delivery store at -20C. Dilute in PBS (pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

## Target


<b>Target Name:</b>	hyaluronan binding protein 2
<b>Alternative Name:</b>	FSAP HABP HGFAL NMTC5 PHBP
<b>Database Link:</b>	<a href="#">NP_004123</a> <a href="#">Entrez Gene 3026 Human</a> <a href="#">Entrez Gene 0 Mouse</a> <a href="#">Entrez Gene 0 Rat</a> <a href="#">Entrez Gene 0 Monkey</a> <a href="#">Entrez Gene 0 Dog</a>
<b>Function:</b>	The protein encoded by this gene is an extracellular serine protease that binds hyaluronic acid and is involved in cell adhesion. The encoded protein is synthesized as a single chain, but then undergoes an autoproteolytic event to form the functional heterodimer. Further autoproteolysis leads to smaller, inactive peptides. This protease is known to cleave urinary plasminogen activator, coagulation factor VII, and the <a href="http://www.amsbio.com">AMS BIO   www.amsbio.com</a>   <a href="mailto:info@amsbio.com">info@amsbio.com</a>

alpha and beta chains of fibrinogen, but not prothrombin, plasminogen, or the gamma chain of fibrinogen. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2010].


## Validation Data



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HABP2 (RC206545, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HABP2.(1:500)



Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-HABP2 monoclonal antibody (1:500)



Immunohistochemical staining of paraffin-embedded Human breast tissue within the normal limits using anti-HABP2 mouse monoclonal antibody.(Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3 min, TA809299)(1:2000)

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

**UK & Rest of the World**  
184 Park Drive, Milton Park  
Abingdon OX14 4SE, UK  
T: +44 (0)1235 828 200  
F: +44 (0) 1235 820 482

**North America**  
1035 Cambridge Street,  
Cambridge, MA 02141  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

**Germany**  
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880

**Switzerland**  
Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
T: +41(0) 91 604 55 22  
F: +41(0) 91 605 17 85