


Product datasheet for TP723918

SARS-CoV-2 Envelope Recombinant Protein with His and Avi tag

Product data:

Product Type:	Recombinant Proteins
Description:	SARS-CoV-2 Envelope Recombinant Protein with His and Avi tag
Species:	SARS-CoV-2
Expression Host:	E. coli
Tag:	N-6xHis/AVI
Purity:	> 95% by SDS-PAGE.
Buffer:	Supplied as a 0.22 µm filtered solution in 20mM Tris,250mM NaCl,0.5%TritonX-100,pH 8.0.
Endotoxin:	< 1.0 EU/µg of the protein by LAL method.
Locus ID:	43740570

Product images:



Recombinant 2019-nCoV envelope Protein with His and Avi tag was determined by SDS-PAGE with Coomassie Blue, showing a band at 12 kDa.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85